



A CLOSER LOOK AT THE

# DIGIRAD X-ACT+

The Digirad X-ACT+ is a solid-state, triple-head, SPECT/FAC camera. With built-in attenuation correction, X-ACT+ offers superior quality.

**SOLID STATE**  
Low-Dose Imaging



**SPECT/FAC**

Fluorescence Attenuation Correction

**OPEN DESIGN**

All Patient Sizes

**ERGONOMIC**

Easy Patient Access

## Next Generation

# ATTENUATION CORRECTION

## 2 in 1

Captures both the Emission and Transmission scans without patient repositioning.



High specificity, PET-like diagnostic value lowers false positives.

## FAC

Low Dose Fluorescence Attenuation Correction is built-in and fully integrated.

SEE THE DIFFERENCE



## UNMATCHED PATIENT EXPERIENCE



Comfortable upright seat



Fast Scan Times



Reduces Motion



Easy-On Easy-Off



No Claustrophobia



All sizes can be accommodated

The X-ACT+ has three solid-state detectors and uses rapid imaging detector geometry.




Experience True **Solid-State**

- ▶ Optimized Image Quality
- ▶ Compact Size
- ▶ TruACQ™ Personalized Imaging

## EASY INSTALL

The X-ACT+ allows for a simple and cost effective installation. The camera requires little space and is essentially plug-and-play.




**Small Size**

X-ACT+ fits in a standard exam room



**No Lead Lining**

X-ACT+ does not require lead lined rooms



**Standard Floor**

The camera weight allows for a standard floor



**110 Power**

No need for special power modifications



Digirad makes healthcare convenient. We are the nationwide leader in delivering diagnostic expertise on an as needed, when needed, where needed basis.

For more information about the Digirad X-ACT+ visit [www.digirad.com](http://www.digirad.com)

[www.digirad.com](http://www.digirad.com) | 800.947.6134 | [info@digirad.com](mailto:info@digirad.com)